**LUNGENE® COVID-19 Antigen Rapid Test** 

REF ICOV5002-B025

For in vitro diagnostic use only.

## [INTENDED USE]

The COVID-19 Antigen Rapid Test is a lateral flow immunoassay intended for the qualitative detection of SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab, nasal swab or oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

Results are for the identification of SARS-CoV-2 nucleocapsid antigen. Antigen is generally detectable in nasopharyngeal swab, nasal swab and oropharyngeal swab during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine clinical status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for clinical diagnosis or to exclude SARS-CoV-2 infection. Clinical judgment and diagnostic results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary for patient management.

The COVID-19 Antigen Rapid Test is intended for medical professionals or trained operators who are proficient in performing lateral flow tests. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation.

## [SUMMARY]

The novel coronavirus (SARS-CoV-2) belongs to the B genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be the main source of infection. Based on the current epidemiological investigation, the incubation period is 2 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

## [PRINCIPLE]

The COVID-19 Antigen Rapid Test is a lateral flow immunoassay based on the principle of lateral flow immunoassay. Sample, test cassette, SARS-CoV-2 nucleocapsid monoclonal antibody conjugate with color microsphere as detector and spray on conjugation pad. During the test, SARS-CoV-2 antigen in the specimen interacts with SARS-CoV-2 antibody conjugate with color microsphere, making antigen-antibody label complex, this complex binds on the membrane via capture antibody. After the reaction, wash the tube with water. Then, transfer SARS-CoV-2 nucleocapsid protein monoclonal antibody. A colored test line (T) would be visible in the result window if SARS-CoV-2 antigen is present in the specimen. At least one T line suggests a negative result. The control line (C) is used for procedural control and should always appear if the test procedure is performed properly.

## [WARNINGS AND PRECAUTIONS]

- For in vitro diagnostic use only.
- For healthcare professionals and individuals trained in point of care settings.
- Do not use this product as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform clinical status of COVID-19.
- Do not use this product after the expiration date.
- Please read all the information in this leaflet before performing the test.
- The test cassette should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test cassette should be discarded according to federal, state and local regulations.

## [COMPOSITION]

**Materials Provided**

- One test cassette, each cassette with desiccant in individual foil pouch
- 25 Extraction Reagent Tubes, tube sealed with foil film containing 0.3 mL of extraction reagent
- 25 Sterile Swabs: single use swab for specimen collection
- 1 Pack Insert

**Materials Required but not Provided**

- Timer

## [STORAGE AND STABILITY]

• Store as packaged in the sealed pouch at temperature (4-30 °C or 40-46 °F). The kit is stable within the expiration date printed on the labeling.

• Once open the pouch, the test should be used within one hour. Prolonged exposure to hot and humid environment will cause product deterioration.




• The LOT and expiration date were printed on the labeling.

## [SPECIMEN]

Specimens obtained early during symptom onset will contain the highest viral titers; specimens obtained later during symptom onset will contain lower viral titers. Results will be compared to an RT-PCR assay. Inadequate specimen collection, improper specimen handling, and/or transport may yield false results; therefore, training in specimen collection is highly recommended due to the importance of specimen quality to accurate test performance.

Prepare the extraction reagent tube according to the Test Procedure and use the sterile swab provided in the kit for specimen collection.

## Nasopharyngeal Swab Specimen Collection


## Specimen Transport and Storage

Do not return the swab to the original swab packaging. Freshly collected specimens should be processed as soon as possible, but no later than one hour after specimen collection. Specimen collected may be stored at 2-8 °C for no more than 24 hours; Store at +4 °C for a long time, but avoid repeated freeze-thaw cycles.

## [TEST PROCEDURE]


Note: Allow the test cassette, reagents and specimens to equilibrate to room temperature (15-30 °C or 59-86 °F) prior to testing.

- Carefully open the foil film on the extraction reagent tube. Do not let the extraction reagent flow out.
- Put the extraction reagent tube on the work station.
- Sampling refers to section Specimen Collection.

- 4. Remove the test cassette from the sealed pouch.
- 5. Reverse the specimen extraction reagent tube, holding the tube upright, transfer 3 drops (approximately 100 µL) slowly to the specimen well (S) of the test cassette, then start the timer.

6. Wait for colored lines to appear. Interpret the test results at 15 minutes. Do not read results after 20 minutes.

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5. Follow Steps 4 - 6 of the Direct Swab Test Procedure above.

## [INTERPRETATION OF RESULTS]

Positive: Two lines appear. One colored line appears at the control region (C), and another colored line appears at the test region (T), regardless of the intensity of the test line.

Negative: One colored line appears at the control region (C), and no line appears at the test region (T).

Control line fails to appear: Insufficient specimen volume or quality. Insufficient specimen volume or quality. The most likely reason for control line failure is that the specimen is too dilute. Repeat the test using a new test cassette.

If the problem persists, discontinue using this test cassette.

All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.

The used test cassette should be discarded according to federal, state and local regulations.

## [QUALITY CONTROL]

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify test performance.

## [LIMITATIONS]

The product is intended to provide a qualitative detection. The intensity of the test line (T) is not intended to correlate to the concentration of the antigen of the specimen.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions.

A physician must interpret the results in conjunction with the patient's history, physical findings, and other diagnostic procedures.

A negative result can occur if the quantity of SARS-CoV-2 antigens present in the specimen is below the detection threshold of the assay, or the virus has undergone minor amino acid mutation(s) in the target epitope region recognized by the monoclonal antibodies utilized in the test.

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